

**Product Name :**  
Cloud and Pour Point Testing Apparatus

**Product Code :**  
TEST5077RFQLAB-0002



**Description :**

Cloud and Pour Point Testing Apparatus for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

For determining the cloud point of an oil provided it is transparent in layer 40 mm in thickness and for the pour point of any petroleum oil.

It consists of a cylindrical glass jar fitted with a cork, a water tight cylindrical jacket with a disc of cork or felt at the bottom, a ring around the glass jar and a bath with cover.

We are a leading manufacturers of Cloud and Pour Point Testing Apparatus testing equipment for a range of industries, including materials science testing lab. We offer full range of Cloud and Pour Point Testing Apparatus testing equipment, including Cloud and Pour Point Testing Apparatus testing kit, Cloud and Pour Point Testing Apparatus testing machine, Cloud and Pour Point Testing Apparatus testing tools, Cloud and Pour Point Testing Apparatus test equipment, Cloud and Pour Point Testing Apparatus test machines, Cloud and Pour Point Testing Apparatus for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Cloud and Pour Point Testing Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

## Testing lab instruments India

