

**Product Name :**  
Vicat Apparatus Modified

**Product Code :**  
TEST5054RFQLAB-0002



#### Description :

Vicat Apparatus Modified for Testing Lab Equipment Tender Suppliers India

#### Technical Specification :

For determination of normal consistency of hydraulic cement, gypsum plaster and lime and setting time of cement, gypsum plaster and gypsum concrete.

The Apparatus consists of the following :-

A frame with a vertically moving bearing rod.

Mould in the form of frustum of a cone 70 mm bottom dia, 60 mm top dia and 40 mm high with base plate.

Needle 1 mm dia x 50 mm long.

Plunger 10 mm dia x 50 mm long.

Plunger 19 mm dia x 44.4 mm long.

Plunger 12.5 mm dia x 25 mm long.

One brass weight.

We are a leading manufacturers of Vicat Apparatus Modified testing equipment for a range of industries, including materials science testing lab. We offer full range of Vicat Apparatus Modified testing equipment, including Vicat Apparatus Modified testing kit, Vicat Apparatus Modified testing machine, Vicat Apparatus Modified testing tools, Vicat Apparatus Modified test equipment, Vicat Apparatus Modified test machines, Vicat Apparatus Modified for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Vicat Apparatus Modified testing equipment that helps customers build better, safer, and more sustainable infrastructure.



## Testing lab instruments India