

**Product Name :**  
Automatic Vicat Apparatus

**Product Code :**  
TEST5054RFQLAB-0003



**Description :**

Automatic Vicat Apparatus

**Technical Specification :**

Suitable for determination of consistency and setting time of standard cement paste, hydraulic cement & Gypsum plaster conforming to IS, EN, ASTM, BS Standards

Fully Automatic test cycle

Electronic sensor to record the data

Storing of test results, settings and data

Large graphic display to follow the test in real time

RS 232 output for connection to PC and printer

Facility to program test cycles as test start delay, dropping rate auto-mode, free or driven drop

Suitable for operation on 230 V AC, 50Hz , Single phase

We are a leading manufacturers of Automatic Vicat Apparatus testing equipment for a range of industries, including materials science testing lab. We offer full range of Automatic Vicat Apparatus testing equipment, including Automatic Vicat Apparatus testing kit, Automatic Vicat Apparatus testing machine, Automatic Vicat Apparatus testing tools, Automatic Vicat Apparatus test equipment, Automatic Vicat Apparatus test machines, Automatic Vicat Apparatus for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Automatic Vicat Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.



**Testing lab instruments India**