

**Product Name :**  
Southard Viscosimeter

**Product Code :**  
TEST5054RFQLAB-0005



**Description :**

Southard Viscosimeter

**Technical Specification :**

For determination of volume yield of quick lime and adjusting the consistency of lime putty mixture.

The equipment consists of the following :-

One 50 mm internal dia metal cylinder mounted on wooden stand and having a close fitting piston with 65 mm working stroke.

A metal bridge for measuring slump.

One lime putty density vessel 63 mm I.D. and 80 mm deep.

We are a leading manufacturers of Southard Viscosimeter testing equipment for a range of industries, including materials science testing lab. We offer full range of Southard Viscosimeter testing equipment, including Southard Viscosimeter testing kit, Southard Viscosimeter testing machine, Southard Viscosimeter testing tools, Southard Viscosimeter test equipment, Southard Viscosimeter test machines, Southard Viscosimeter for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Southard Viscosimeter testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

## Testing lab instruments India

