

Product Name :
Jolting Apparatus

Product Code :
TEST5060RFQLAB-0008



Description :

Jolting Apparatus

Technical Specification :

Used for molding of three standard specimens of 40 mm x 40 mm x 160 mm size from Portland and Pozzolana cement mortar for determination of their transverse strengths. The equipment consists of the following: –
A motorised Jolting Machine, as per IS 10078 specifications, which gives 60 jolts per minute to the rectangular table through a reduction gear and cam drive. A pre-set type counter with a micro switch stops the machine automatically after 60 jolts. The rectangular table has guide pins for properly locating and clamping the moulds.
A three compartment mould for preparing three specimens of 40 mm x 40 mm x 160 mm size.
A three compartment hopper which can be mounted on the mould for pouring the mortar in to the mould.
Suitable for operation on 220 V, 50 Hz, single phase supply.

We are a leading manufacturers of Jolting Apparatus testing equipment for a range of industries, including materials science testing lab. We offer full range of Jolting Apparatus testing equipment, including Jolting Apparatus testing kit, Jolting Apparatus testing machine, Jolting Apparatus testing tools, Jolting Apparatus test equipment, Jolting Apparatus test machines, Jolting Apparatus for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Jolting Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India