

**Product Name :**  
Clamped Bosses Proving Rings

**Product Code :**  
TEST5022RFQLAB-0001



**Description :**

Clamped Bosses Proving Rings for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

These Proving rings have clamped bosses and are available in the range of 0.5kN – 100kN. repeatability and accuracy of all the rings

Capacity

0.5 kN

1.0 kN

2 kN

2.5 kN

5 kN

10 kN

20 kN

25 kN

50 kN

100 kN

**Integral Proving Rings**

These Proving rings have integral loading bosses forged and machined along with the ring. These bosses can not shift their position to cause mis-alignment which may result in loss of accuracy These rings are manufactured up to 3000kN capacity.

Capacity

---

0.5 kN  
1.0 kN  
2 kN  
2.5 kN  
5 kN  
10 kN  
20 kN  
25 kN  
50 kN  
100 kN  
150 kN  
200 kN  
250 kN  
500 kN  
1000 kN  
2000 kN  
3000 kN

We are a leading manufacturers of Clamped Bosses Proving Rings testing equipment for a range of industries, including materials science testing lab. We offer full range of Clamped Bosses Proving Rings testing equipment, including Clamped Bosses Proving Rings testing kit, Clamped Bosses Proving Rings testing machine, Clamped Bosses Proving Rings testing tools, Clamped Bosses Proving Rings test equipment, Clamped Bosses Proving Rings test machines, Clamped Bosses Proving Rings for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Clamped Bosses Proving Rings testing equipment that helps customers build better, safer, and more sustainable infrastructure.



**Testing lab instruments India**