

**Product Name :**  
Plastic Limit Apparatus

**Product Code :**  
TEST5001RFQLAB-0002



**Description :**

Plastic Limit Apparatus for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Conforming (Part V) The apparatus consists of :  
Glass plate 200 mm x 150 mm x 3 mm thick  
Porcelain evaporating dish 120 mm dia.  
Brass rod 3 mm dia, 100 mm long.  
Flexible spatula with 80 mm long, 20 mm wide blade.  
Moisture Cans 50 mm dia, 12 Nos.

We are a leading manufacturers of Plastic Limit Apparatus testing equipment for a range of industries, including materials science testing lab. We offer full range of Plastic Limit Apparatus testing equipment, including Plastic Limit Apparatus testing kit, Plastic Limit Apparatus testing machine, Plastic Limit Apparatus testing tools, Plastic Limit Apparatus test equipment, Plastic Limit Apparatus test machines, Plastic Limit Apparatus testing Instruments. We are suppliers of materials testing tools equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Plastic Limit Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

## Testing lab instruments India

