

**Product Name :**  
Marshall Stability Test Apparatus

**Product Code :**  
TEST5122RFQLAB-0001



#### Description :

Marshall Stability Test Apparatus for Testing Lab Equipment Tender Suppliers India

#### Technical Specification :

Max load: 30KN  
Speed: 50±5mm/min  
Proving ring: 30KN  
Micrometer: 0-10mm Accuracy?0.01mm  
Power supply: AC220V 50Hz  
Motor capacity: 550W

We are a leading manufacturers of Marshall Stability Test Apparatus testing equipment for a range of industries, including materials science testing lab. We offer full range of Marshall Stability Test Apparatus testing equipment, including Marshall Stability Test Apparatus testing kit, Marshall Stability Test Apparatus testing machine, Marshall Stability Test Apparatus testing tools, Marshall Stability Test Apparatus test equipment, Marshall Stability Test Apparatus test machines, Marshall Stability Test Apparatus for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Marshall Stability Test Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

## Testing lab instruments India

