

Product Name :
Automatic Wheel-Tracker Immersion Test

Product Code :
TEST5125RFQLAB-0002



Description :

Automatic Wheel-Tracker Immersion Test for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Wheel tracker is intended for evaluating the anti-rutting ability of asphalt mixture at high temperature and helping with the ratio design of asphalt mixture.

Working mode: Water immersion test and non-immersion test

Temperature measurement accuracy: $\pm 0.1^{\circ}\text{C}$

Ambient temperature control accuracy: $\pm 1^{\circ}\text{C}$

Number of temperature measurement channels: 2

Displacement measurement accuracy: $\pm 0.01\text{mm}$

Displacement measurement range: 0~30mm

Temperature setting range: 40°C ~ 80°C

We are a leading manufacturers of Automatic Wheel-Tracker Immersion Test testing equipment for a range of industries, including materials science testing lab. We offer full range of Automatic Wheel-Tracker Immersion Test testing equipment, including Automatic Wheel-Tracker Immersion Test testing kit, Automatic Wheel-Tracker Immersion Test testing machine, Automatic Wheel-Tracker Immersion Test testing tools, Automatic Wheel-Tracker Immersion Test test equipment, Automatic Wheel-Tracker Immersion Test test machines, Automatic Wheel-Tracker Immersion Test for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Automatic Wheel-Tracker Immersion Test testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India