


**Product Name :**  
Polishing Grinding Machine

**Product Code :**  
TEST5158RFQLAB-0002



**Description :**

Polishing Grinding Machine for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Grinding and polishing machine can effectively replace previously manual grinding and polishing processes, which greatly saves labor, improves processing quality, and elevates sample preparation efficiency.

Grinding disc diameter: 230mm  
Abrasive paper diameter: 230mm  
Rotating speed: 900r/min  
Electromotor: 0.2kW, 380V, 50Hz

We are a leading manufacturers of Polishing Grinding testing equipment for a range of industries, including materials science testing lab. We offer full range of Polishing Grinding testing equipment, including Polishing Grinding testing kit, Polishing Grinding testing machine, Polishing Grinding testing tools, Polishing Grinding test equipment, Polishing Grinding test machines, Polishing Grinding for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Polishing Grinding testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

## Testing lab instruments India

