

Product Name :
Automatic Water Distiller

Product Code :
TEST5176RTLBT-0019



Description :

Automatic Water Distiller

Technical Specification :

Automatic Water Distiller

All-in-one, automatically-controlled distillation water distiller made of stainless steel with an inbuilt storage tank.

Specifications:

Electronic control system (automatic water and power cut-off when tank is full and automatic start, Stainless steel heaters)

Thermostatic low water cut off system

Original design, which allows to produce distilled water with temp. about 35-40° C without separate chiller

Electronic leakage protection system with solenoid valve

Easy-access evaporator, condenser-tank for effortless maintenance (cleaning, etc)

Easily replaceable condenser

Automatic level control in the evaporator

Visual control of the distillate level in the tank (through glass level-tube on the front panel)

All parts that contact water or steam are made of high-quality stainless steel

Stainless steel tubular heater with extended resource

Thermostatic low-water cut-off safety system.

Output: 4 liters/ hour

Water consumption better: 70 liters/ hour or

Storage tank Capacity tank: 8 liters
Distillate conductivity: 2.5 μ S/m
Power Supply: 230V 50-60Hz, single phase



Testing lab instruments India