

**Product Name :**  
Steel Transverse Flexural Accessory For Steel Tester

**Product Code :**  
TEST5176RTLBT-0042



**Description :**

Steel Transverse Flexural Accessory For Steel Tester

**Technical Specification :**

Steel Transverse Flexural Accessory For Steel Tester  
 Type Of Testing Machines: Flexure Testing Machine  
 Material: Mild Steel  
 Finishing: Polished  
 Size: 4X6feet  
 Automation Grade: Automatic  
 Frequency: 50 Hz  
 Surface Finish: Polished  
 Maximum Load (KN): 200  
 Maximum distance between bearer (mm): 800  
 Bearer length (mm): 120  
 Bearer diameter (mm): 50.

