

**Product Name :**  
Orbital Shaker Digital and Horizontal

**Product Code :**  
TEST5178RFQLAB-0003



**Description :**

Orbital Shaker Digital and Horizontal

**Technical Specification :**

Orbital Shaker, Digital, and Horizontal

With digital speed display, interchangeable shaker platform

Speed range (0-300 rpm), accept shaking loads up to 15kg.

Two different shaking functions: back and forth, and orbital

Adjustable shaking speed and range

Timer function for flexibility

Corrosion resistant stainless steel housing,

Easy to clean and maintain

Large mounting surface and big load capacity of up to 15kg

Ideal for vessels with more than 250 ml capacity, like round flasks, Erlenmeyer flasks, culture flasks, or culture bottles

Widely used for bacterial cultures, fermentation, hybridization, biochemical reactions, enzyme or tissue research

Shaking motion can be set in up to 6 different ways.

---

## Testing lab instruments India

