

Product Name :
Cyclo Mixer

Product Code :
TEST5178RFQLAB-0006



Description :

Cyclo Mixer

Technical Specification :

Cyclo Mixer

For delicate oliquid mixing with bi- directional switch, speed regulation, and interchangeable mixing used for variety of tubes.

Steel construction, 210-240V, 50/60Hz AC power supply

With LCD display.

Capacity 250 mL.

Made from top quality components in sync with top-notch technology.

Features: Efficient performance, longer service life.

Variable speed mixer to eliminate hand mixing. Speed regulator controls the degree of vibration. Suitable for touch/continuous operations.

Supplied with Interchangeable adaptors for different applications.

Designed for mixing liquids in Schools, Laboratories & Factories

Touch/ Continuous Operation mode Selection through bi-directional Switch

Speed Regulation through knob provided on the control panel

Interchangeable mixing heads for use with variety of tubes

Supplied with all interchangeable mixing heads.



Testing lab instruments India