

**Product Name :**  
Vacuum Oven Vacuum Drying Cabinet

**Product Code :**  
TEST5178RFQLAB-0007



**Description :**

Vacuum Oven Vacuum Drying Cabinet

**Technical Specification :**

Vacuum Oven (Vacuum Drying Cabinet)

Electric thermostatic air drying oven is produced by advanced laser and numerical control processing equipment. It is used for drying and sterilization.

Capacity: 70L.

Temperature Range: RT + 10 ~ 250 °, the temperature fluctuation degree, plus or minus 1.0 °, the temperature accuracy: plus or minus 0.1 °

Temperature controller: digital computer P.I.D temperature controller, with timing function, accurate and reliable temperature control;

Alarm function: independent temperature limit alarm system, automatic interruption when the temperature exceeds the limit, no accident occurs in the safe operation of the experiment (optional);

Box material: the liner material is mirror stainless steel shell material is quality steel plate spraying plastics;

Timing range: 1-9999min;

Working environment: + 5 ~ 40 °;

Power source: AC 220V±10% 50HZ.

---

## Testing lab instruments India

