


Product Name :
Manual Capsule Filler (Size 0)

Product Code :
TEST5178RFQLAB-0010



Description :

Manual Capsule Filler (Size 0)

Technical Specification :

Manual Capsule Filler (Size 0)

Size 0 high precision 20 holes manual capsule filler/encapsulation machine, suitable for the separated capsule. Raw material is imported acrylic which is healthy, good transparency, durable, low weight.

The screws & pins are made of stainless steel and meet the hygiene requirement GMP standards.

The cap sheet adopt one-layer design and polished angles to prevent bleeding.

The middle sheet adopt special flute design of high precision, which can avoid the breakage of capsules when joining them.

The body sheet adopt the specially-dealt screws and springs, which can make no noise when pressing downwards. The detachable and undetachable body sheet is optional.

The encapsulation sheet adopt special design, which can increase the speed of capsules collection.

The tamping tool is design to one-layer ,angles and pins being polished.

The powder guard with a U-design powder guard to avoid loss of powders.

