

**Product Name :**  
Disintegration Tester

**Product Code :**  
TEST5178RFQLAB-0016



**Description :**

Disintegration Tester

**Technical Specification :**

Disintegration tester: is available with up to 4 racks for baskets with either 3 or 6 test stations.  $55 \pm 2$  mm stroke height, Accuracy  $\pm 0.2^\circ$  C, temperature range  $30-50^\circ$  C. Up to 100 samples per test.

Disintegration Testers, the test stations (one to four) are driven individually. In order to guarantee outstanding quality and the highest ease of use each test station is controlled by its own keypad.

After programming the required test run time, the basket rack assembly with its loaded samples is automatically lowered into the media. On completion of the preset run time or when the stop button is pressed, the basket rack is raised from the media to prevent further disintegration of the samples (optional).

The temperature displayed on the LED display is the one measured by the standard external temperature sensor, either located inside the water bath or the test media. The optional type "B" basket rack assembly with three larger test tubes can be fitted to all disintegration testers.

