

**Product Name :**  
Vibrating Table

**Product Code :**  
TEST5180RFQLAB-0003



**Description :**

Vibrating Table

**Technical Specification :**

Vibrating Table

For compacting concrete specimen in the laboratory vibrating table mounted on a steel stand, supplied with clamp assembly.

Proper compaction of cement concrete while casting specimens for compression testing is essential to achieve higher compressive strength. Vibrating Table is ideally suited for this purpose.

The table top is suitable to hold cube moulds and has stops along its edges to prevent moulds from sliding off the table during operation. The specially designed vi bro motor for operating the vibrator.

Table size: 60g x 800 mm

Amplitude: 0.3-0.6 mm

Vibrating frequency: 3000 vibrations/ min

Power: 220 V, &.3 kW

---