

**Product Name :**  
Laboratory Mixer (Mortar Mixer)

**Product Code :**  
TEST5180RFQLAB-0005



**Description :**

Laboratory Mixer (Mortar Mixer)

**Technical Specification :**

Laboratory Mixer (Mortar Mixer)

This mixer is designed to mix mortars and cement pastes to the requirements of the above standard. The mixing paddle has a planetary motion and is driven by a motor with a microprocessor-based speed and program controller.

The mixer can be operated either in an automatic or manual mode. When the mixer is used in the manual mode, the two mixing speeds can be changed by means of a rocker switch, without switching off the motor. In the automatic mode any one of the pre-set mixing programmes may be selected.

Rotation speed of mixing blade revolution(r/min):

four: 140±5 62±5

high: 285±10 125±10

Mixing bowl capacity: 5 L

Width of mixing blade: 135 mm

Power Supplied: AC 220V, 300KW

---

## Testing lab instruments India

