

**Product Name :**  
Aggregate Impact Value

**Product Code :**  
TEST5194RFQLAB-0037



**Description :**

Aggregate Impact Value

**Technical Specification :**

Aggregate impact value  
Aggregate Impact Value Test Apparatus is used for determining the aggregate impact value. The Sturdy Construction consists of a base and support columns to form a rigid framework around the quick release trigger mechanism to ensure an effective free fall of the hammer during test.  
The free fall can be adjusted through  $380 \pm 5$ mm. The hammer is provided with a locking arrangement.  
Aggregate impact value test apparatus is supplied complete with a cylindrical measure of 75mm dia x 50mm depth, an automatic blow counter and a tamping rod.  
Size: LxWxH: 500x300x980mm